# *Polypektom* **ENDOX**

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
## ENDOx – the partially insulated snare for precise removal without causing thermal damage to the organ walls

The speed at which polyps or lesions are removed en bloc is of great consequence for the successful closure of vessels. A slower incision results in a more effective outcome; however, if the incision is made too slowly, thermal damage to the organ wall may occur.

The ENDOx polypectomy snare takes a new approach. One half of the symmetrically opening snare is electrically insulated and HF-surgically inactive. The first cut and the cutting phase are therefore extra effective. The low use of HF electricity minimises both thermal damage to the muscular layer and serous membrane as well as thermal artefacts.

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### Endoscopic check of the cut speed thanks to the markings on the resection snare


If the cut is made too fast it may cause bleeding, while if the cut is too slow it may lead to thermal blood vessel damage. To minimise these risks the insulated part of the ENDOx snare is equipped with markings, like a guidewire, which can be seen by the endoscope thanks to the transparent catheter.

This allows the endoscope user and assistants to see clearly on the monitor whether or how fast the snare closes and whether it cuts.

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
When using conventional polypectomy snares on large polyps (> 1cm), long first cut delays can cause big and deep blood vessel damage (to the coagulation plus devitalisation zones) (see arrow).

## Less thermal damage to the muscular layer thanks to the shorter first cut delay

The HF-surgically active cutting wire of the ENDOx snare is around 50% shorter than on similarly-sized, non-insulated snares. The first cut delay is also much shorter as a result. Furthermore, the risk of cutting towards the organ wall with the ENDOx snare is only half of the risk compared to conventional polypectomy snares. Consequently the ENDOx snare is particularly suitable for removing large polyps.

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When using ENDOx resection snares the first cut delay is relatively short even with large polyps (> 1cm) and as a result the thermal blood vessel damage (to the coagulation plus the adjacent devitalisation zones) is only relatively small and not deep (see arrow).

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## Ergonomic 3-finger-handle for controlled cutting

For the precise application of the snare and cutting, the ENDOx snare is equipped with an ergonomically-shaped handle. Axial and rotational movements are transferred very precisely from the handle to the snare. At the same time, when the snare touches the organ you can feel it clearly with your hand. The catheter is also separated from the handle, allowing the snare to be precisely opened and closed.

- No friction or jamming
- Minimum mechanical clearance
- HF connection at the side
- The high-quality multi-contact plug avoids loose connections

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## Technical Data

|                                                                           |                                            |
|---------------------------------------------------------------------------|--------------------------------------------|
| Opening and closing of the effector                                       | symmetrical                                |
| Rotating the snare                                                        | at the catheter axis                       |
| Maximum opening width of snare                                            | ca. 3- 4 cm                                |
| Diameter of the cutting wire                                              | 0,38 mm                                    |
| Length of the effective wire                                              | 50% of similarly-sized conventional snares |
| Length of the catheter                                                    | 200 cm                                     |
| External diameter of the catheter                                         | 2,3 mm                                     |
| Plugs available on ENDOx for cables to<br>connect on a HF surgical device | 4 mm, male                                 |

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Feinwerktechnik GmbH